

Lys-to-Lys Protein-Protein Conjugation Kit (Small Scale 100-500 μ g)

Product Number: CCR-1008

Description

Lys-to-Lys Protein-Protein Conjugation Kit (Cat. # 1008) provides sufficient reagents to perform two protein-protein conjugation reactions. Start to finish, any two lysine containing proteins (100-500 μ g) in a volume of 100 μ L (1-5 mg/mL) can be efficiently conjugated in less than 3 hours.

Kit Contents:

Component	Amount
Tetrazine-PEG ₅ -NHS ester(Tz)	2 x 0.5 mg
<i>Trans</i> -Cyclooctene-PEG ₄ -NHS ester (TCO)	2 x 0.5 mg
DMSO	5 mL
BupH TM Saline Buffer Pack	1 pack
Zeba TM Spin Columns (0.5 mL)	8

Introduction

Lys-to-Lys Protein-Protein Conjugation kit provides all the necessary reagents to perform two protein-protein conjugation reactions. Conjugates are formed by an inverse-electron demand Diels-Alder cycloaddition reaction between *trans*-cyclooctene (TCO) and tetrazine (Tz) functional groups, a bioorthogonal click reaction characterized by exceptional kinetics ($k > 800 \text{ M}^{-1}\text{s}^{-1}$) and selectivity (Figure 1). The Tz/TCO reaction pair represents the most powerful tool available for catalyst-free protein-protein bioconjugation. Fast kinetics combined with exquisite specificity enables rapid conjugation (30-60 min) of lysine containing proteins to each other at low concentrations (e.g. 5-10 μ M) with > 99% conversion of the limiting protein to conjugate in mild buffered media (e.g. PBS pH 7.5). Other features include long-term reactive stability of TCO and Tz functional groups on modified proteins stored in aqueous buffered media (e.g. maintaining > 90% reactivity after 1 month at 4 $^{\circ}$ C, pH 7.5). This stability avoids time-sensitive timing of reactions associated with classical chemistries based on thiol/maleimide chemistry, providing worry-free conjugation results.

Important Information

- NHS esters are moisture-sensitive. Avoid moisture condensation by allowing product to come to room temperature before opening. Prepare working stock solutions immediately before use.
- For NHS ester reactions avoid buffers containing primary amines (e.g. Tris, glycine).




Figure 1. Schematic representation of protein-protein conjugation chemistry.

Protein Requirements

- Protocol requires 100-500 μg of each protein in a volume of 0.1 mL (1-5 mg/mL)
- Proteins must be highly purified and their molecular weight known
- Proteins must have available primary amines (e.g. lysine residues)

Additional Information

- When conjugating proteins at low concentrations (e.g. $<5 \mu\text{M}$), we recommend using a 4-5 fold molar excess over the limiting protein.
- Whenever possible, maintain the protein used in excess at higher concentration than limiting protein. For example, if IgG is the limiting protein at 1 mg/mL keep the excess protein at 4-5 mg/mL for best results.
- Final protein-protein conjugates can be purified by a number of methods including size exclusion, ion exchange, affinity, or hydrophobic interaction chromatography. For certain applications, protein-protein conjugates can be used without further purification.

Materials Required but Not Provided

- UV-VIS spectrophotometer
- Microcentrifuge capable of handling 1.5 mL tubes
- Quartz semi-micro cuvette (50-100 μL)
- 1.5 mL microfuge tubes
- Pipettes and tips (P-10, P-100, P-1000)
- Ultrapure water (e.g. 18 M Ω -cm)
- Beaker, stir bar, and 6 N NaOH
- pH meter

Material Preparation

A. BupH™ Buffer Preparation

1. Dissolve BupH™ dry-blend buffer pack (provided) into 480 mL ultrapure water. Adjust the pH of the solution to 7.5 ± 0.05 with 6N NaOH. Adjust the final volume to 500 mL with ultrapure water. For long-term storage sterile-filter the solution. Do not add sodium azide or Proclin 300 preservatives as these reagents interfere with protein determination (A280).

B. Protein Preparation

1. If protein #1 (100-500 μg) is lyophilized and free of exogenous amines (e.g. glycine or Tris), resuspend in 100 μL BupH buffer (pH 7.5) to obtain a 1-5 mg/mL solution. Proceed to Tetrazine Labeling of Protein #1 as described in Section E.
2. If protein #2 (100-500 μg) is lyophilized and free of exogenous amines (e.g. glycine or Tris), resuspend in 100 μL BupH buffer (pH 7.5) to obtain a 1-5 mg/mL solution. Proceed to TCO Labeling of Protein #2 as described in Section F.
3. If either or both purified proteins (100-500 μg) are already in 100 μL buffer (e.g. PBS), for optimal results buffer exchange into BupH (pH 7.5) using the spin columns provided prior to Tz or TCO labeling. See Step C and Step D for buffer exchange of proteins.

C. Spin Column Equilibration into BupH (pH 7.5) prior to Buffer Exchange

1. Twist off the column's bottom closure and loosen the cap. Place each 0.5 mL spin column into a Clean 1.5 mL microfuge tube.
2. Centrifuge column at 1,500 x g for 2 minutes to remove storage solution. Place a pen mark on the side of the column where the compacted resin is slanted upward. Place column in centrifuge with the mark facing away from the center of the rotor in all subsequent centrifugation steps.
Note-resin will appear white in color and compacted after centrifugation.
3. Add 0.3 mL BupH buffer (pH 7.5) to the top of each spin column, replace the cap and loosen.
4. Centrifuge at 1,500 x g for 2 minutes to remove buffer.
5. Repeat steps 3 and 4 two additional times, discarding buffer from collection tube after each spin.
6. Transfer equilibrated spin column (resin appears white and dry) into a new 1.5 mL microfuge tube and immediately proceed to buffer exchange of protein.

D. Buffer Exchange of Protein

1. Buffer exchange protein into BupH (pH 7.5) equilibrated spin column by slowly applying 100 μL protein solution to the center of equilibrated resin bed.
2. Centrifuge at 1,500 x g for 2 minutes. Retain the eluate at bottom of 1.5 mL collection tube.
3. Protein is now buffer exchanged.

Tetrazine Labeling of Protein #1 (Limiting Protein)

1. Determine volume DMSO required to dissolve Tetrazine-PEG₅-NHS reagent provided (0.5 mg) by referring to the calculations in Appendix A.

2. Add required volume DMSO to tetrazine reagent, vortex for 2 minutes to completely dissolve.
3. Add 5 μL of tetrazine (20-fold molar excess) to 100 μL buffer exchanged protein #1 (limiting protein).
4. Allow reaction to proceed for 60 minutes at room temperature.
5. Remove excess reagent from labeled protein using a new buffer exchange spin column equilibrated in BupH (pH 7.5) by following Steps C and D. Retain final Tz-labeled protein #1 at bottom of collection tube.
6. Determine concentration of Tz-modified protein #1 (mg/mL) by measuring A280. Refer to calculations in Appendix A. **Note**-remove a 10 μL aliquot of Tz-modified protein and dilute into 90 μL BupH H (e.g. 1:10), then measure A280 using a semi-micro quartz cuvette (50-100 μL). Alternately, a Bradford or BCA protein assay can be performed using a 10 μL aliquot of the labeled protein.
7. Tetrazine-labeled protein is now ready for conjugation.

TCO Labeling of Protein #2 (Excess Protein)

1. Determine volume DMSO required to dissolve TCO-PEG₄-NHS reagent provided (0.5 mg). Refer to Appendix B for calculations.
2. Add required volume DMSO to TCO reagent, vortex for 2 minutes to completely dissolve.
3. Add 5 μL of TCO (20-fold molar excess) to 100 μL buffer exchanged protein #2.
4. Allow reaction to proceed for 60 minutes at room temperature.
5. Remove excess reagent from labeled protein using a new buffer exchange spin column equilibrated in BupH (pH 7.5) by following Steps C and D. Retain final TCO-labeled protein #2 at bottom of collection tube.
6. Determine concentration of TCO-labeled protein #2 (mg/mL) by measuring A280. Refer to calculations in Appendix B. **Note**- remove a 5 μL aliquot of TCO-modified protein and dilute into 95 μL BupH H (e.g. 1:20), then measure A280 using a semi-micro quartz cuvette (50-100 μL). Refer to calculations in Appendix B. Alternately, a Bradford or BCA protein assay can be used to determine protein concentration if the protein's E1% or molar extinction coefficient are unknown.
7. TCO-labeled protein is now ready for conjugation.

Protein-Protein Conjugation

1. Select the desired protein-protein stoichiometry for your conjugation reaction (e.g. 1:3)
Note-typical conditions use a 1.2-5 fold molar excess of the excess protein over the limiting protein.
2. Using the selected stoichiometry, calculate volume of TCO-modified protein (excess protein) to add to Tz-modified protein (limiting protein). Refer to calculations in Appendix C.
3. Mix required volume of TCO-modified protein with Tz-modified protein.
4. Allow the conjugation reaction to proceed for 60 minutes at room temperature.
5. Store protein-protein conjugates at 4°C until ready for use or purification.

Troubleshooting

Problem	Possible Cause	Solution
No conjugation of tetrazine and TCO modified proteins	One or more proteins are not properly labeled with TCO or Tetrazine	Confirm purity and concentration of proteins prior to labeling process. Buffer exchange proteins into BupH (pH 7.5) if necessary.
	NHS-ester hydrolyzed	Allow product to equilibrate to room temperature before opening.
		Avoid buffers that contain primary amines such as Tris and glycine. Buffer exchange proteins before labeling if necessary.
Excess labeling reagent improperly desalted	Remove excess un-reacted tetrazine and TCO reagents by desalting.	
Low conjugation of TCO and Tetrazine labeled proteins	Suboptimal reaction conditions	Optimize conjugation conditions by altering molar excess
		Confirm proper concentration of protein #1-Tz and protein #2-TCO prior to conjugation (i.e. 1-5 mg/mL)