


Methyltetrazine-Amine, HCl salt

Product Number: CCR-1011

Product Name: Methyltetrazine-Amine, HCl salt

Chemical Structure:



CAS Number: 1345955-28-3

Chemical Composition: C₁₀H₁₁N₅

Molecular Weight: 201.23 (237.83 HCl salt)

Solubility: DMSO, DMF

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Slightly pink solid

Shelf Life: At least 12 months at -20°C

For Research Use Only