


Methyltetrazine-PEG4-Amine, HCl salt

Product Number: CCR-1012

Product Name: Methyltetrazine-PEG4-Amine, HCl salt

Chemical Structure:



CAS Number: 1802908-05-9

Chemical Composition: C₁₇H₂₅N₅O₄

Molecular Weight: 363.41 (HCl salt)

Solubility: DMSO, DMF, DCM, THF, Chloroform

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Pink solid

Shelf Life: At least 12 months at -20°C

For Research Use Only