


Methyltetrazine-PEG4-Azide

Product Number: CCR-1014

Product Name: Methyltetrazine-PEG4-Azide

Chemical Structure:



CAS Number: N/A

Chemical Composition: C₁₇H₂₃N₇O₄

Molecular Weight: 389.40

Solubility: DMSO, DMF, DCM, MeOH

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Red solid

Shelf Life: At least 12 months at -20°C

For Research Use Only