


Cy5 Methyltetrazine

Product Number:	CCR-1019
Product Name:	Cy5 Methyltetrazine
CAS Number:	1801695-58-8
Molecular Weight:	934.11
Solubility:	Water, DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C.
Appearance:	Blue solid
Shelf Life:	At least 12 months at -20°C


Chemical Structure:



Optical Properties

λ_{abs}	647 nm
λ_{abs}	663 nm
ϵ_{max}	255,000 M ⁻¹ cm ⁻¹
Φ_{fl}	0.21
CF ₂₆₀	0.05
CF ₂₈₀	0.05
Flow cytometry laser line	633 or 635 nm
Microscopy laser line	633 or 635 nm
Spectrally Similar Dyes	Alexa Fluor 647 CF Dye 647, DyLight 649

Absorption and Emission Spectra



For Research Use Only