


Biotin-PEG4-Methyltetrazine

Product Number: CCR-1024

Product Name: Biotin-PEG4-Methyltetrazine

Chemical Structure:



Molecular Weight: 589.71

Chemical Composition: $C_{27}H_{39}N_7O_6S$

CAS Number: N/A

Solubility: DMSO, DMF, MeOH, DCM, THF

Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.

Appearance: Red solid or oil.

Shelf Life: At least 12 months at $-20^{\circ}C$.

For Research Use Only