


## Biotin-PEG3-Maleimide

Product Number: CCR-1029  
Product Name: Biotin-PEG3-Maleimide  
Chemical Structure:



CAS Number: 1431618-70-0  
Molecular Weight: 597.72  
Chemical Composition:  $C_{27}H_{43}N_5O_8S$   
Solubility: DMSO, DMF, MeOH  
Storage: Upon receipt store at  $-20^{\circ}C$ . Product shipped at ambient temperature.  
Appearance: Orange solid  
Shelf Life: At least 12 months at  $-20^{\circ}C$

**For Research Use Only**