


## Cy5.5 DBCO

Product Number:	CCR-1046
Product Name:	Cy5.5 DBCO
CAS Number:	N/A
Molecular Weight:	1175.36
Solubility:	Water, DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C.
Appearance:	Blue solid
Shelf Life:	At least 12 months at -20°C


### Chemical Structure:



### Optical Properties

$\lambda_{abs}$	678 nm
$\lambda_{abs}$	694 nm
$\epsilon_{max}$	190,000 M <sup>-1</sup> cm <sup>-1</sup>
$\Phi_{fl}$	0.28
CF <sub>260</sub>	0.10
CF <sub>280</sub>	0.13
Spectrally Similar Dyes	Alexa Fluor 680, IRDye 680RD CF Dye 680, DyLight 680

### Absorption and Emission Spectra



For Research Use Only