


## TBTA

Product Number: CCR-1061  
Product Name: TBTA  
Alternative Name: Tris[(1-benzyl-1*H*-1,2,3-triazol-4-yl)methyl]amine  
Chemical Structure:



CAS Number: 510758-28-8  
Molecular Weight: 530.62  
Chemical Composition: C<sub>30</sub>H<sub>30</sub>N<sub>10</sub>  
Solubility: DMSO, DMF, DCM, THF  
Storage: Upon receipt store at 4°C. Product shipped at ambient temperature.  
Appearance: Off-white to grey solid  
Shelf Life: At least 12 months at 4°C

**For Research Use Only**