


Methyltetrazine-PEG4-Maleimide

Product Number: CCR-1068

Product Name: Methyltetrazine-PEG4-Maleimide

Chemical Structure:



CAS Number: N/A

Molecular Weight: 514.53

Chemical Composition: $C_{24}H_{30}N_6O_7$

Solubility: DMSO, DMF, THF, DCM, MeOH

Storage: Upon receipt store at -20°C . Product shipped at ambient temperature.

Appearance: Red solid

Shelf Life: At least 12 months at -20°C

For Research Use Only