


## Methyltetrazine-PEG4-NHS Ester

Product Number: CCR-1069

Product Name: Methyltetrazine-PEG4-NHS Ester

Chemical Structure:



CAS Number: 1802907-92-1

Chemical Composition:  $C_{24}H_{31}N_4O_9$

Molecular Weight: 533.53

Solubility: DMSO, DMF, MeOH, DCM, THF

Storage: Upon receipt store at  $-20^{\circ}\text{C}$ . Product shipped at ambient temperature.

Appearance: Red solid

Shelf Life: At least 12 months at  $-20^{\circ}\text{C}$

**For Research Use Only**