


Methyltetrazine-Propylamine, HCl Salt

Product Number: CCR-1071

Product Name: Methyltetrazine-Propylamine, HCl Salt

Chemical Structure:



CAS Number: N/A

Molecular Weight: 281.74 (HCl salt)

Chemical Composition: C₁₂H₁₆ClN₅O (HCl salt)

Solubility: DMSO, DMF, THF, DCM

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Slightly pink solid

Shelf Life: At least 12 months at -20°C

For Research Use Only