


TAMRA Acid

Product Number:	CCR-1079
Product Name:	TAMRA Acid
Alternative Names:	5(6)-Carboxytetramethylrhodamine
Chemical Structure:	



CAS Number:	98181-63-6
Molecular Weight:	430.45
Abs/Em:	553/575 nm
E _{max} :	92,000cm ⁻¹ M ⁻¹
Chemical Composition:	C ₂₅ H ₂₂ N ₂ O ₅
Solubility:	DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance:	Dark Red solid
Shelf Life:	At least 12 months at -20°C

For Research Use Only