


## TCO-PEG4-Acid

Product Number: CCR-1088

Product Name: TCO-PEG4-Acid

Chemical Structure:



CAS Number: N/A

Chemical Composition:  $C_{20}H_{35}NO_8$

Molecular Weight: 417.49

Solubility: DMSO, DMF, DCM, EtOAc, MeOH

Storage: Upon receipt store at  $-20^{\circ}\text{C}$ . Product shipped at ambient temperature.

Appearance: Colourless to slightly yellow oil

Shelf Life: At least 12 months at  $-20^{\circ}\text{C}$

**For Research Use Only**