


Cy5-TCO

Product Number: CCR-1089

Product Name: Cy5-TCO

Chemical Structure:



CAS Number: N/A

Chemical Composition: C₄₅H₆₀N₄O₁₂S₃

Molecular Weight: 945.17

Abs/Em: 647/633nm

E_{max}: 251,000cm⁻¹M⁻¹

Solubility: Water, DMSO, DMF, MeOH

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Blue solid

Shelf Life: At least 12 months at -20°C

For Research Use Only