


## Bis-Sulfone NHS Ester

Product Number: CCR-1104

Product Name: Bis-Sulfone NHS Ester

Chemical Structure:



CAS Number: 899452-51-8

Molecular Weight: 597.66

Chemical Composition: C<sub>29</sub>H<sub>27</sub>NO<sub>9</sub>S<sub>2</sub>

Solubility: DMSO, DMF, THF, DCM

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: White amorphous solid

Shelf Life: At least 12 months at -20°C

**For Research Use Only**