


## Desthiobiotin-PEG3-Azide

Product Number: CCR-1107

Product Name: Desthiobiotin-PEG3-Azide

Chemical Structure:



CAS Number: 1306615-47-3

Molecular Weight: 414.50

Chemical Composition: C<sub>18</sub>H<sub>34</sub>N<sub>6</sub>O<sub>5</sub>

Solubility: DMSO, DMF, MeOH

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Colourless to slightly yellow oil

Shelf Life: At least 12 months at -20°C

**For Research Use Only**