


DBCO-PEG4-Desthiobiotin

Product Number: CCR-1108

Product Name: DBCO-PEG4-Desthiobiotin

Chemical Structure:



CAS Number: N/A

Molecular Weight: 719.87

Chemical Composition: $C_{39}H_{61}N_5O_8$

Solubility: DMSO, DMF, DCM, THF, MeOH

Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.

Appearance: Colourless to slightly yellow oil

Shelf Life: At least 12 months at $-20^{\circ}C$

For Research Use Only