


PC Biotin-DBCO

Product Number: CCR-1120

Product Name: PC Biotin-DBCO

Chemical Structure:



CAS Number: N/A

Molecular Weight: 1002.14

Chemical Composition: $C_{50}H_{61}N_4O_{13}S$

Solubility: DMSO, DMF, DCM, THF, MeOH

Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.

Appearance: Slightly yellow amorphous solid

Shelf Life: At least 12 months at $-20^{\circ}C$

For Research Use Only