


## Methyltetrazine-Acid

Product Number: CCR-1125

Product Name: Methyltetrazine-Acid

Chemical Structure:



CAS Number: N/A

Molecular Weight: 230.22

Chemical Composition:  $C_{11}H_{10}N_4O_2$

Solubility: DMSO, DMF, THF, MeOH

Storage: Upon receipt store at  $-20^{\circ}C$ . Product shipped at ambient temperature.

Appearance: Pink solid

Shelf Life: At least 12 months at  $-20^{\circ}C$

**For Research Use Only**