


Methyltetrazine-NHS Ester

Product Number: CCR-1128

Product Name: Methyltetrazine-NHS Ester

Chemical Structure:



CAS Number: 1644644-96-1

Molecular Weight: 327.29

Chemical Composition: $C_{40}H_{48}N_2O_{14}S_2$

Solubility: DMSO, DMF, THF, DCM, Chloroform

Storage: Upon receipt store at -20°C . Product shipped at ambient temperature.

Appearance: Pink solid

Shelf Life: At least 12 months at -20°C

For Research Use Only