


Dde Biotin-Azide

Product Number: CCR-1136

Product Name: Dde Biotin-Azide

Chemical Structure:



CAS Number: N/A

Molecular Weight: 695.37

MW Left Behind: 100.07

Chemical Composition: C₃₂H₅₃N₇O₈S

Solubility: DMSO, DMF, MeOH

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Glass-like solid

Shelf Life: At least 12 months at -20°C

For Research Use Only