


## Dde Biotin-DBCO

Product Number: CCR-1138

Product Name: Dde Biotin-DBCO

Chemical Structure:



CAS Number: N/A

Molecular Weight: 872.08

MW left behind: 276.32

Chemical Composition:  $C_{32}H_{50}N_4O_8S$

Solubility: DMSO, DMF,  $CH_2Cl_2$ , THF

Storage: Upon receipt store at  $-20^\circ C$ . Product shipped at ambient temperature.

Appearance: Glass-like solid

Shelf Life: At least 12 months at  $-20^\circ C$

**For Research Use Only**