


## Tetrazine-PEG<sub>5</sub>-NHS Ester

Product Number:	CCR-1143
Product Name:	Tetrazine-PEG <sub>5</sub> -NHS Ester
Alternative Name:	H-Tetrazine-PEG5-NHS Ester
Chemical Structure:	



CAS Number:	1682653-80-0
Molecular Weight:	604.61
Chemical Composition:	C <sub>27</sub> H <sub>36</sub> N <sub>9</sub> O <sub>10</sub>
Solubility:	DMSO, DMF, MeOH, DCM
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance:	Red amorphous solid to red oil
Shelf Life:	At least 12 months at -20°C

**For Research Use Only**