


5-Ethynyl-2'-Deoxyuridine

Product Number: CCR-1149
Product Name: 5-Ethynyl-2'-Deoxyuridine
Alternative Name: EdU

Chemical Structure:



CAS Number: 61135-33-9
Molecular Weight: 252.23
Chemical Composition: $C_{11}H_{12}N_2O_5$
Solubility: Water, DMSO, DMF
Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.
Appearance: White crystalline
Shelf Life: At least 12 months at $-20^{\circ}C$

For Research Use Only