


N-(4-pentynoyl)-glucosamine tetraacylated

Product Number: CCR-1155
Product Name: N-(4-pentynoyl)-glucosamine tetraacylated (Ac₄GlcNAI)
Chemical Structure:



CAS Number: 1361993-37-4
Chemical Composition: C₁₉H₂₅NO₁₀
Molecular Weight: 427.40
Solubility: DMSO, DMF, ME OH, EtOAc
Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance: Off-white to grey amorphous solid
Shelf Life: At least 12 months at -20°C

For Research Use Only