


## PC DBCO-NHS Ester

Product Number: CCR-1160

Product Name: PC DBCO-NHS Ester

Chemical Structure:



CAS Number: N/A

Chemical Composition:  $C_{47}H_{55}N_5O_{16}$

Molecular Weight: 945.96

Solubility: DMSO, DMF, MeOH, DCM, THF

Storage: Upon receipt store at  $-20^{\circ}\text{C}$ . Product shipped at ambient temperature.

Appearance: Yellow amorphous solid

Shelf Life: At least 12 months at  $-20^{\circ}\text{C}$

**For Research Use Only**