


PC Methyltetrazine-NHS Ester

Product Number: CCR-1163

Product Name: PC Methyltetrazine-NHS Ester

Chemical Structure:



CAS Number: N/A

Molecular Weight: 785.75

Chemical Composition: $C_{35}H_{43}N_7O_{14}$

Solubility: DMSO, DMF, MeOH, DCM, THF

Storage: Upon receipt store at -20°C . Product shipped at ambient temperature.

Appearance: Yellow Oil

Shelf Life: At least 12 months at -20°C

For Research Use Only