


Biotin Picolyl Azide

Product Number: CCR-1167

Product Name: Biotin Picolyl Azide

Chemical Structure:



CAS Number: N/A

Molecular Weight: 622.74

Chemical Composition: $C_{27}H_{42}N_8O_7S$

Solubility: DMSO, DMF, MeOH

Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.

Appearance: Off-white to slightly orange amorphous solid

Shelf Life: At least 12 months at $-20^{\circ}C$

For Research Use Only