


FAM Picolyl Azide

Product Number:	CCR-1180
Product Name:	FAM Picolyl Azide
Molecular Weight:	580.55
Solubility:	DMSO, DMF
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance:	Yellow to orange solid
Shelf Life:	At least 12 months at -20°C



Optical Properties:

λ_{abs}	490 nm
λ_{abs}	510 nm
ϵ_{max}	76,000 M ⁻¹ cm ⁻¹
Φ_{fl}	0.24
CF ₂₆₀	0.24
CF ₂₈₀	0.26
Flow cytometry laser line	488 nm
Microscopy laser line	488 nm
Spectrally Similar Dyes	Alexa Fluor 488 CF Dye 488, DyLight 488

For Research Use Only