


## Cy7 Picolyl Azide

Product Number:	CCR-1183
Product Name:	Cy7 Picolyl Azide
Molecular Weight:	1175.34
Solubility:	Water, DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance:	Dark green solid
Shelf Life:	At least 12 months at -20°C



<b>Optical Properties:</b>	$\lambda_{abs}$	753 nm
	$\lambda_{abs}$	775 nm
	$\epsilon_{max}$	255,000 M <sup>-1</sup> cm <sup>-1</sup>
	$\Phi_{fl}$	0.16
	CF <sub>260</sub>	0.05
	CF <sub>280</sub>	0.05
	Spectrally Similar Dyes	Alexa Fluor 750 IR Dye 750, DyLight 750

**For Research Use Only**