


Tetrazine-Sulfo-NHS Ester

Product Number: CCR-1191

Product Name: Tetrazine-Sulfo-NHS Ester

Chemical Structure:



CAS Number: N/A

Molecular Weight: 413.31

Chemical Composition: C₁₄H₁₀N₅O₇S

Solubility: DMSO, DMF, Water

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Pink solid

Shelf Life: At least 12 months at -20°C

For Research Use Only