


Methyltetrazine-Sulpho-NHS Ester

Product Number: CCR-1193
Product Name: Methyltetrazine-Sulpho-NHS Ester
Chemical Structure:



CAS Number: N/A
Molecular Weight: 429.34
Chemical Composition: C₁₅H₁₂N₅NaO₇S
Solubility: Water, DMSO, DMF
Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance: Pink solid
Shelf Life: At least 12 months at -20°C

For Research Use Only