


## TAMRA Tetrazine

Product Number:	CCR-1196
Product Name:	TAMRA Tetrazine
Molecular Weight:	846.92
Solubility:	Water, DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance:	Red solid
Shelf Life:	At least 12 months at -20°C
Optical Properties:	



$\lambda_{abs}$	553nm
$\lambda_{abs}$	565nm
$\epsilon_{max}$	86,000 M <sup>-1</sup> cm <sup>-1</sup>
$\phi_{fl}$	0.09
Flow Cytometry Laser Line	532, 555 or 568nm
Microscopy Laser Line	532 or 555nm
Spectrally Similar Dyes	Alexa Fluor ®555, CF Dye® 555, DyLight ®549

**For Research Use Only**