


Cy3 Tetrazine

Product Number:	CCR-1204
Product Name:	Cy3 Tetrazine
Molecular Weight:	1331.36
Solubility:	Water, DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance:	Blue solid
Shelf Life:	At least 12 months at -20°C



Optical Properties:

λ_{abs}	555nm
λ_{abs}	565nm
ϵ_{max}	155,000 M ⁻¹ cm ⁻¹
ϕ_{fl}	0.09
CF ₂₆₀	0.05
CF ₂₈₀	0.05
Flow Cytometry Laser Line	532, 555 or 568nm
Microscopy Laser Line	532 or 555nm
Spectrally Similar Dyes	Alexa Fluor ®555, CF Dye® 555, DyLight ®549

For Research Use Only