


## PC Biotin-PEG3-NHS Ester

Product Number: CCR-1225  
Product Name: PC Biotin-PEG<sub>3</sub>-NHS Ester  
Chemical Structure:



CAS Number: N/A  
Molecular Weight: 840.89  
Chemical Composition: C<sub>36</sub>H<sub>52</sub>N<sub>6</sub>O<sub>15</sub>S  
Solubility: DMSO, DMF, DCM, THF  
Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.  
Appearance: Slightly yellow amorphous solid to oil  
Shelf Life: At least 12 months at -20°C

**For Research Use Only**