


DBCO-C6-NHS Ester

Product Number: CCR-A102

Product Name: DBCO-C6-NHS Ester

Chemical Structure:



CAS Number: 1384870-47-6

Molecular Weight: 430.45

Chemical Composition: $C_{25}H_{22}N_2O_5$

Solubility: DMSO, DMF, DCM, THF, Chloroform

Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.

Appearance: Grey solid

Shelf Life: At least 6 months at $-20^{\circ}C$

For Research Use Only