


DBCO-PEG5-NHS Ester

Product Number: CCR-A102P
Product Name: DBCO-PEG5-NHS Ester
Chemical Structure:



CAS Number: 1378531-80-6
Molecular Weight: 693.74
Chemical Composition: $C_{18}H_{16}N_2O$
Solubility: DMSO, DMF, DCM, THF, $CHCl_3$
Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.
Appearance: Yellow to slightly orange oil
Shelf Life: At least 6 months at $-20^{\circ}C$

For Research Use Only