


DBCO-PEG4-Amine

Product Number: CCR-A103P

Product Name: DBCO-PEG4-Amine

Chemical Structure:



CAS Number: 1255942-08-5

Molecular Weight: 523.62

Chemical Composition: C₂₉H₃₇N₃O₆

Solubility: DMSO, DMF, THF, CH₂Cl₂, CHCl₃

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Yellow Oil

Shelf Life: At least 6 months at -20°C

For Research Use Only