


DBCO-PEG4-Biotin

Product Number: CCR-A105

Product Name: DBCO-PEG4-Biotin

Chemical Structure:



CAS Number: 1255942-07-4

Molecular Weight: 749.92

Chemical Composition: C₃₉H₅₁N₅O₈S

Solubility: DMSO, DMF, THF, DCM, MeOH

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: White amorphous solid or oil

Shelf Life: At least 12 months at -20°C

For Research Use Only