


DBCO-PEG4-Maleimide

Product Number: CCR-A108P

Product Name: DBCO-PEG4-Maleimide

Chemical Structure:



CAS Number: 1480516-75-3

Molecular Weight: 674.74

Chemical Composition: C₃₆H₄₂N₄O₉

Solubility: DMSO, DMF, CH₂Cl₂, CHCl₃, THF

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Yellow Oil

Shelf Life: At least 6 months at -20°C

For Research Use Only