


Disulphide DBCO Biotin

Product Number:	CCR-A112
Product Name:	Disulphide DBCO Biotin
Alternative Name:	DBCO-S-S-PEG3-Biotin
Chemical Structure:	



CAS Number:	N/A
Molecular Weight:	869.12
Chemical Composition:	$C_{42}H_{56}N_6O_3S_3$
Solubility:	DMSO, DMF
Storage:	Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.
Appearance:	Glass-like amorphous solid
Shelf Life:	At least 12 months at $-20^{\circ}C$

For Research Use Only