


DBCO-Sulpho-NHS Ester

Product Number: CCR-A124
Product Name: DBCO-Sulpho-NHS Ester

Chemical Structure:



CAS No.: 1400191-52-7
Molecular Weight: 532.50
Chemical Composition: C₂₅H₂₁N₂NaO₈S
Solubility: Water, DMSO, DMF
Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance: White to grey solid
Shelf Life: At least 12 months at -20°C

For Research Use Only