


Carboxyrhodamine 110-DBCO

Product Number:	CCR-A127
Product Name:	Carboxyrhodamine 110-DBCO
Alternative Name:	ADIBO-PEG ₄ -Carboxyrhodamine 110, DBCO-PEG ₄ -Carboxyrhodamine 110

Chemical Structure:



Molecular Weight:	880.96
Abs/Em:	501/526 nm
E _{max} :	74,000 M ⁻¹ cm ⁻¹
Chemical Composition:	C ₅₀ H ₅₀ N _n O ₁₀
Solubility:	DMSO, DMF
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Shelf Life:	At least 12 months at -20°C

For Research Use Only