


DBCO-PEG4-DBCO

Product Number: CCR-A128

Product Name: DBCO-PEG4-DBCO

Chemical Structure:



CAS Number: N/A

Molecular Weight: 854.92

Chemical Composition: $C_{50}H_{54}N_4O_9$

Solubility: DMSO, DMF, DCM, THF, CH_2Cl_2 , $CHCl_3$

Storage: Upon receipt store at $-20^{\circ}C$. Product shipped at ambient temperature.

Appearance: Grey solid

Shelf Life: At least 12 months at $-20^{\circ}C$

For Research Use Only