


TAMRA-DBCO

Product Number: CCR-A131

Product Name: TAMRA-DBCO

Alternative Name: ADIBO-PEG₄-TAMRA, DBCO-PEG₄-TAMRA

Chemical Structure:



Molecular Weight: 937.07

Abs/Em: 555/571 nm

E_{max}: 92,000cm⁻¹M⁻¹

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Solubility: DMSO, DMF

Appearance: Red solid

Shelf Life: At least 12 months at -20°C

For Research Use Only