


## DBCO-NHS Ester

Product Number: CCR-A133

Product Name: DBCO-NHS Ester

Chemical Structure:



CAS Number: 1353016-71-3

Molecular Weight: 402.40

Chemical Composition:  $C_{23}H_{18}N_2O_5$

Solubility: DMSO, DMF, DCM, THF, Chloroform

Storage: Upon receipt store at  $-20^{\circ}C$ . Product shipped at ambient temperature.

Appearance: White solid

Shelf Life: At least 12 months at  $-20^{\circ}C$

**For Research Use Only**