


TCO-PEG4-NHS Ester

Product Number: CCR-A137

Product Name: TCO-PEG4-NHS Ester

Chemical Structure:



CAS Number: 1621096-79-4

Molecular Weight: 514.57

Chemical Composition: $C_{24}H_{38}N_2O_{10}$

Solubility: DMSO, DMF, DCM, THF, Chloroform

Storage: Upon receipt store at -20°C . Product shipped at ambient temperature.

Appearance: Yellow Oil

Shelf Life: At least 12 months at -20°C

For Research Use Only