


## Cy3 DBCO

Product Number:	CCR-A140
Product Name:	Cy3 DBCO
CAS Number:	N/A
Molecular Weight:	983.18. 92
Solubility:	Water, DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C.
Appearance:	Red solid
Shelf Life:	At least 12 months at -20°C


Chemical Structure:



### Optical Properties

$\lambda_{\text{abs}}$	550 nm
$\lambda_{\text{em}}$	570 nm
$\epsilon_{\text{max}}$	150,000 M <sup>-1</sup> cm <sup>-1</sup>
$\Phi_{\text{fl}}$	0.10
CF <sub>260</sub>	0.05
CF <sub>280</sub>	0.05
Flow cytometry laser line	532 or 555 nm
Microscopy laser line	532 or 555 nm
Spectrally Similar Dyes	Alexa Fluor 555 CF Dye 555, DyLight 555

Absorption and Emission Spectra



For Research Use Only