


Carboxyrhodamine 110-Azide

Product Number:	CCR-AZ105
Product Name:	Carboxyrhodamine 110-Azide
Alternative Name:	Azide-PEG3-5(6)-Carboxyrhodamine 110, Azide-PEG3- Rhodamine Green, Azide Fluor 488

Chemical Structure:



CAS Number:	N/A
Molecular Weight:	575.59
Chemical Composition:	C ₂₉ H ₃₁ N ₆ O ₇
Abs/Em:	501/525 nm
E _{max} :	74,000cm ⁻¹ M ⁻¹
Appearance:	Brick red solid
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Shelf Life:	At least 12 months at -20°C

For Research Use Only